Optima

Optima

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## General information

Please read the instructions for use carefully each time before you use the products. The products must only be used by doctors, individuals instructed in their proper use and individuals entrusted with patient care.

**CAUTION!** For use on a single patient only!

### 1. Intended use

**OPTIMA** tracheostomy tubes are used to hold open a tracheostoma and permit the patient to breathe. Some versions allow phonation and, as a result, speech.

### 2. Product description

#### 2.1. General

**OPTIMA** tracheostomy tubes are produced from a soft and elastic, biocompatible plastic. The flexibility of the material ensures good manoeuvrability combined with high wearing comfort. The eyelets on the sides of the neck flange are used to retain the inserted tracheostomy tube in place with the help of a neck strap.

The 15-mm universal connector is permanently fixed to the flange of the outer tube. The inner cannula clicks into the 15-mm connector. Inflatable versions also have a particularly thin-walled low-pressure cuff for sealing the tube flush against the wall of the trachea. A control balloon with a non-return valve is used to inflate/deflate the low-pressure cuff and to monitor the applied cuff pressure.

**OPTIMA** tracheostomy tubes are always supplied with an obturator with an atraumatic tip and 2 inner cannulas, which can be easily removed if and as required, e.g., for cleaning or in case of respiratory distress.

The fenestrated inner cannula in the tube versions with a **speaking option (Voice)** is red in colour. Since the inner cannula is visible on the side away from the patient even when it has been inserted, the user and/or nursing personnel are alerted immediately to the fact that a fenestrated inner cannula is in use.

The standardised 15-mm universal connector ensures the secure connection of accessories with a 15-mm connection such as an HME, a speaking valve or ventilation machines to all **OPTIMA** tracheostomy tubes.

All of the materials used in their manufacture are biocompatible and comply with the applicable legal standards.

#### 2.2. Product versions

up-to-date overview of the tracheostomy tubes and accessories currently available can be found in our product catalogue, which you can request from us free of charge or view online at [www.primed-halberstadt.de](http://www.primed-halberstadt.de).

A table giving the dimensions of the **OPTIMA** tracheostomy tubes can be found at the end of the instructions for use.

### 2.3. Cannula pass

The cannula pass is enclosed with every package containing a tracheostomy tube. After sticking the provided label on it, the pass contains information on:

- *Article number*
- *Length*
- *Size*
- *Batch no.*

## 3. Period of use

We recommend exchanging the tube:

- after 29 days

(from the first use respectively, including all periods in which the tube is temporarily not in use).

We strongly recommend always having at least two replacement tubes on hand to ensure uninterrupted treatment at all times.

## 4. Indications

For ensuring an optimal fit and air supply, the corresponding **OPTIMA** tracheostomy tube must always be selected by the attending doctor.

### General

- Keeping open of airways in case of functional or mechanical obstruction of the upper respiratory tract
- Long-term ventilation

### OPTIMA versions without cuff (Basic, Voice):

- For spontaneously breathing patients who rely on a tube.

### OPTIMA versions with cuff (Basic Cuff, Voice Cuff, Suction, Suction Voice):

- For patients in whom an access to the airways via a tracheostoma with sealing off of the trachea is necessary.

### OPTIMA versions VOICE, with sieving/fenestration (Voice, Voice Cuff, Suction Voice):

- For tracheotomised patients with a fully or partially preserved larynx. These versions enable the patients to speak.

### OPTIMA versions SUCTION, with suction function (Suction, Suction Voice):

- For patients in whom suctioning from the subglottic region is indicated.

## 5. Contraindications

- Do not insert **OPTIMA** tracheostomy tubes in patients with abnormal anatomy or pathology of the upper airways as there is a risk of partial or total obstruction of the airways.
- Do not use **OPTIMA** tracheostomy tubes with a speaking valve/decannulation cap in laryngectomised (voice box surgically removed) patients (**risk of suffocation!**), **Exception: Voice prosthesis/shunt valve.**
- Due to the increased risk of aspiration, **VOICE** versions must never be used in patients with recurring aspiration pneumonia.
- Use in paediatrics.

- Trachea located unusually far from the surface of the body (e.g., in obese patients), possibility of using longer tubes.

## 6. Side effects

Pressure wounds, stenoses, necroses, skin irritations, granulation tissue, tussive irritation (urge to cough), dysphagia and bleeding.

## 7. Warnings /Safety precautions

1. Do not exert any force on the tracheostomy tube (**risk of damage**).
2. When the patient is moved, he/she must not lie on the control balloon.
3. The filling tube must not lie between the inner and outer tube when the inner cannula is replaced (**risk of damage**).
4. Only remove items from the 15-mm universal connector using approved tools (e.g., disconnecting wedge).
5. Only insert sieved and fenestrated tubes under medical supervision during ventilation.
6. The cuff should not come into contact with aerosols or ointments containing lidocaine.
7. Thin cuffs display a certain degree of moisture permeability and, as a result, condensation may collect in the cuff (minor quantities are no cause for concern).
8. Tracheostomy tubes with **SUCTION function** may only be used in patients with an increased risk of bleeding (e.g., when taking anticoagulants) if the attending doctor deems it appropriate due to the higher risks during suction procedures.
9. The **VOICE versions** can be connected to a speaking valve equipped with a corresponding connection via the integrated, standardised 15-mm universal connector of the outer tube (**OBSERVE: The manufacturer's instructions for use!**). The red one of the two inner cannulas supplied has an oval opening (fenestration) extending over the sieving of the outer tube so that air can flow freely for speaking when the patient exhales.
10. **VOICE versions:** These tube versions are only recommended for patients with heavy secretions or a tendency to granulation tissue and clogging if they are regularly checked by a doctor and the intervals between replacement are shorter (weekly, as a rule) because the sieving in the outer tube can intensify the formation of granulation tissue.
11. Only the manufacturer is permitted to shorten lengths, add sieving holes or otherwise modify the tubes, as patients may otherwise be put at risk. Any work on tracheostomy tubes which is not performed properly can lead to serious injury! All liability claims are voided if modifications are carried out by unauthorised individuals.
12. Ensure the tube and the ventilation hose are securely positioned in ventilated patients as dislocation or disconnection can pose a risk of lethal injury!
13. There is a decannulation cap included with the **VOICE versions** for use in preparation for possible decannulation following a temporary tracheostomy. This temporarily stops the air flow through the tube so that the patient becomes accustomed to breathing through his/her mouth/nose again. Decannulation and thus interruption

of the breathing via the inserted tube must only be performed in patients with a preserved larynx and under medical supervision (**always deflate the cuff in cuff versions**). **RISK OF SUFFOCATION!**

14. Do not attempt to reposition a tracheostomy tube with inflated low-pressure cuff while it is being worn! **RISK OF INJURY! DEFLATE CUFF FIRST! Suction off secretions if necessary!**
15. Maintain a safe distance to the tracheostomy tube when using lasers or radiation therapies. **RISK OF FIRE!**
16. Speaking valves must not be inserted when the patient is sleeping.
17. The cuff pressure can change as a result of altitude changes (e.g., in aeroplanes), when nitrous oxide (laughing gas) is used in anaesthesia and when a hand-held pressure gauge is connected or disconnected.
18. Excessive cuff pressure poses the risk of permanent damage to the trachea.
19. Insufficient cuff pressure poses the risk of aspiration.
20. There is a funnel connector included with the **SUCTION versions**. This must only be connected to the suction line. Connections to the cuff pressure line can result in immediate deflation.
21. Irritations, tussive irritation and bleeding are possible when the tubes are inserted and removed.
22. **MRI safety information:**

Optima		
Basic	✓	
Basic Voice	✓	
Cuff		✓
Voice Cuff		✓
Suction		✓
Suction Voice		✓

#### **MR Safe**

There are no restrictions for the use of these tracheostomy tubes during MRI scans. It must be noted that accessories containing metal (e.g., silver speaking valves) must be removed from the tracheostomy tubes prior to the scans.

#### **MR Conditional**

Non-clinical studies have shown that these tracheostomy tubes are only suitable for MRI scans under certain conditions. Patients with these products can undergo an MRI scan safely under the following conditions.

- Accessories containing metal (e.g., silver speaking valves) must be removed from the tracheostomy tubes prior to the scans.
- Static magnetic field of up to 3.0 tesla (T).
- Maximum spatial field gradient of 11 T/m.

- Body coils with quadrature transmission only.
- Secure the tracheostomy tube with a neck strap in order to avoid possible dislocation during the MRI scans.
- In tracheostomy tubes with cuffs, the pilot balloon must be secured in place (e.g. with medical tape) at a safe distance from the region being scanned.

The quality of the MRI results may be affected under certain conditions if the region to be scanned is close to the tracheostomy tube or pilot balloon.

Tracheostomy tubes with cuffs: In non-clinical studies, the visual artefact caused by the pilot balloon extended over up to 100 mm.

## 8. Application

### 8.1. Insertion of the tracheostomy tube

#### 8.1.1. Preparation

The entire tracheostomy tube must be inspected for external damage and loose parts so as to exclude the possibility of injuries and guarantee the necessary sealing. Do not use the tube if you identify anything abnormal.

#### OPTIMA versions with cuff (Basic Cuff, Voice Cuff, Suction, Suction Voice):

The airtightness of the tube and low-pressure cuff should be checked before every insertion. To do so, fill the low-pressure cuff with a pressure of 40 mmHg ( $\approx$  54.4 cmH<sub>2</sub>O or 53.3 hPa). A cuff pressure measuring device is recommended for the test (**OBSERVE: The manufacturer's instructions for use!**).

**If the above-mentioned signs appear, indicating that the tube is not airtight, use of this tube must not be continued under any circumstances. The proper functioning of this tube is no longer assured!**

#### 8.1.2. Insertion

**Tubes with low-pressure cuff must be completely deflated before insertion!**

We recommend sliding a tube compress over the tube shaft to protect the tracheostoma and keep the patient's skin and clothing dry.

We recommend applying stoma oil to the outer tube to reduce the friction on the tracheostomy tube and simplify its insertion into the trachea. This also reduces tussive irritation.

To insert the tracheostomy tubes, first remove the inner cannula and then insert the obturator. Tilt the head backwards gently. Grasp the tube by the flange, using the other hand to spread open the tracheostoma slightly if necessary, then carefully insert the end of the tube into the stoma opening during the inspiration phase (when the patient is breathing in). As you advance the tube, the head is returned to its normal position at the same time and the tube inserted vertically down into the trachea following the arc of the product.

Following insertion, remove the obturator immediately and insert the inner cannula. To fix the inner cannula in position, hold the flange of the outer tube in position with your fingertips and insert the inner cannula into the fixed 15-mm universal connector until it clicks into place audibly. If the patient is inserting the tube independently, he/she should use a mirror to simplify the insertion of the tracheostomy tube. This allows him/her to observe every action. This tube should be fixed in place with the neck strap for secure positioning.

### **8.1.3. Inflating the low-pressure cuff**

Once the tube has been inserted and fixed in place, a cuff pressure meter is used to fill the low-pressure cuff with air. The inflation pressure must be adjusted for the specific requirements and is determined by the treating physician.

Unless otherwise ordered by the doctor, we recommend a cuff pressure of at least 16 to 18 mmHg /22 to 24 cmH<sub>2</sub>O /21 to 24 hPa. If the low-pressure cuff is airtight, there will not be any spontaneous, significant loss of pressure. The pressure should be checked every two hours.

### **8.2. Removing the tracheostomy tube**

**Disconnect the tracheostomy tube from the ventilation system and remove any accessories as well as the neck strap.**

#### **Tracheostomy tubes with suction function (OPTIMA Suction, OPTIMA Suction Voice):**

Always check the cuff pressure and adjust if necessary before performing suction!

When using the **SUCTION versions**, suction the secretion off using the suction system integrated in the tube before removing. Connect the funnel connector to the tube and then switch on the suction device. (**OBSERVE: The manufacturer's instructions for use!**)

**The negative pressure during suction must not exceed 0.2 bar (-200 hPa).** Suction can also be performed via the Luer connector using a syringe. In order to avoid risks and depending on the circumstances of the specific illness, the procedure must be discussed and determined together with the attending doctor. A procedure which is as optimal and gentle as possible must be chosen.

The design of the **SUCTION versions** also allows performance of subglottic flushing (e.g., with saline solution, NaCl 0.9%, exact dosage according to doctor's instructions). This can reduce mucous membrane irritations and soften incrustations. During the flushing procedure, the usual suction process is reversed and the solution is introduced into the trachea above the low-pressure cuff via the integrated suction system. Upon conclusion of the procedure, the solution must be removed again using a syringe or suction device as described above.

#### **OPTIMA versions with cuff, without suction function (Basic Cuff, Voice Cuff):**

In the versions with cuff and without suction function, secretions should be suctioned off first after the inner cannula is removed. This should be done by inserting a suction catheter/line through the tube down to above the tip of the tube with the tracheostomy tube in place. Any secretions which are present can be suctioned off by simultaneously applying suction and deflating the low-pressure cuff. This prevents any aspiration of the secretions.

**⚠ The versions with cuff may not be removed from the tracheostoma until the low-pressure cuff has been completely deflated: Failure to do so can cause injuries. If it is not possible to deflate the cuff and the filling line is not kinked, cut the filling line between the neck flange and the control balloon.**

The patient's head should be tilted back slightly when the tube is removed. To remove the tube, hold it by the sides of the flange. When doing so, proceed extremely carefully so as to avoid causing injuries.

**⚠ In case of unstable tracheostomata or in emergency cases (puncture tracheostoma, dilation tracheostoma), a new tube should be kept ready as the tracheostoma can collapse following removal of the tube and impede the patient's air supply.**

- ⚠** The tightness of the seal and the adequacy of the cuff pressure must be checked before the flushing procedure is carried out. Suction off any secretions above the low-pressure cuff before beginning the flushing procedure.

### 8.3. Exchanging the inner cannulas

Grasp the inner cannula by the recesses on the universal connector and pull it out carefully. Then insert a replacement inner cannula into the fixed 15-mm universal connector of the outer tube until it is heard to click into place audibly. Clean or dispose of the old inner cannula.

- ⚠** Do not twist the inner cannula along its axis when removing it as this may damage it. Do not use force when removing the inner cannula. If the inner cannula cannot be released, remove the entire tube.

## 9. Cleaning/ Disinfection/ Care

### 9.1. Cleaning

Tracheostoma tubes and their accessories must be cleaned regularly, with the cleaning interval always adapted to the individual patient's requirements. Any accessories (speaking valve, HMEs, etc.) must be removed and cleaned separately before the cleaning. The inner cannula is to be removed from the outer cannula.

Use lukewarm water and a mild, pH-neutral washing lotion for the cleaning.

- ⚠** Observe and follow the manufacturer's instructions for use when using cleaning agents.

#### Method

The use of a cannula soaking container with a tray insert simplifies handling (**observe and follow the manufacturer's instructions for use!**). The outer and inner cannulas are immersed in the prepared cleaning solution alongside each other. In versions with a cuff, the cuff must be inflated before cleaning and the control balloon left outside the cleaning bath. Ensure that outer and inner cannulas are completely immersed and there is no air trapped inside. Always place only the outer and inner tubes of one tracheostomy tube in the soaking container to avoid mix-ups. Once the exposure time has elapsed, the outer and inner cannulas must be rinsed thoroughly multiple times under running water.

If there are any stubborn and tough secretion residues which cannot be removed with the cleaning bath, manual cleaning using a special tube cleaning brush is possible. Always insert the cleaning brush into the tube from the tip. The tubes must subsequently be rinsed thoroughly again under running water.

After washing the tube, dry it thoroughly with a clean, lint-free cloth.

There must be no cleaning agent residues remaining on the tube when it is inserted into the tracheostoma.

- ⚠** OPTIMA tracheostomy tubes should only be cleaned with agents approved by the manufacturer.
- ⚠** Do not use any aggressive household cleaners, high-proof alcohol, hydrogen peroxide or products for cleaning dentures! This would pose an acute risk to health. In addition, the tracheostomy tube could be damaged or destroyed.
- ⚠** Never use a dishwasher, steam cooker or microwave oven to clean the product!

⚠ **The product must not be allowed to boil/overheat during the cleaning.**

### **9.2. Disinfection**

Disinfection should only be performed if the attending doctor has ordered it because of the specific clinical picture or if it is indicated by the particular care situation.

#### **Method**

For cold disinfection, the outer and inner cannulas need to be cleaned as described above and then immersed in a disinfectant solution suitable for tracheostomy tubes (e.g., **PRIMEDICLEAN®**) (**Observe and follow the manufacturer's instructions for use!**). The disinfection process is identical to the cleaning process.

There must be no disinfectant residues remaining on the tube when it is inserted into the tracheostoma.

If necessary, obturators can be cleaned and disinfected in the same way as the cannulas.

Before reinserting the tube, check that it is free from damage and functioning perfectly!

Any liquid penetrated into the cuff is to be withdrawn via the control balloon with the help of a syringe.

⚠ **Failure to comply with the cleaning/disinfection specifications will void all guarantee claims.**

### **9.3. Care**

After cleaning, disinfection if applicable, and drying of the tube, the outer surface of the outer tubes, with the exception of versions with low-pressure cuff, should be lubricated by the application of stoma oil or lubricant gel (f.ex. **Optilube®**). This makes them easier to insert. Regular care preserves the malleability of the outer tube's material.

## **10. Storage/ Warehousing**

Tubes which have not yet been used should be stored in a dry environment, protected from direct sunlight and/or heat.

## **11. Disposal**

The devices must be disposed of in accordance with the respectively applicable national regulations governing waste materials.

## **12. Reporting**

All serious events and incidents arising in connection with this medical device must be reported to the manufacturer and the responsible authorities.

## **13. Guarantee**

We guarantee that the delivered object is free from defects in accordance with the scope of our general terms and conditions (GTCs). You can find a copy of our GTCs among our business documentation and on our website at [www.primed-halberstadt.de](http://www.primed-halberstadt.de), etc.

Repairs and other work on the **OPTIMA** tracheostomy tubes may only be carried out by **Primed Halberstadt Medizintechnik GmbH**.

Failure to comply with these terms will void the warranty immediately!

Products are subject to manufacturer's modifications at any time.

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