Diarflex® Pro

- DE Gebrauchsanweisung
- EN Instructions for use
- FR Consignes d'utilisation
- ES Instrucciones de uso
- PT Manual de instruções
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- DA Brugsanvisning
- FI Käyttöohje
- LT Naudojimo instrukcija
- CS Návod k použití
- HU Használati utasítás
- SK Návod na použitie
- BG Инструкции за употреба
- EL Οδηγίες χρήσης
- NO Bruksanvisning
- TR Kullanım kılavuzu
- AR مَدْخِطَسَالَا لِيَالِد

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General information

Read the directions for use carefully each time before using the products. Use of the product is reserved for physicians, trainees and persons entrusted with care.

Suitable for single patient use only.

1. Intended purpose

The **Diarflex® Pro** stool drainage system is intended for draining liquid and semi-liquid stool and guarantees that it will be collected safely for bedridden or immobile and incontinent patients (including irrigation and endoluminal medication application on doctor's orders, as well as taking samples).

2. Product description

The contents and description of the components can be found on the fold-out pages.

3. Useful life

Intended for use for up to 29 consecutive days.

4. Indications

- Complete or partial loss of control over defecation (e.g. bedridden, immobile or incontinent patients)
- Liquid or semi-liquid stools
- Endoluminal drug delivery (in accordance with a doctor's prescription)
- Prevention of contamination of wounds with faeces

5. Contraindications

- Do not use for patients with a known sensitivity or allergy to the materials used in this product
- Do not use if the patient's rectum cannot hold the fill volume of the low pressure balloon or if the anal canal is constricted (anal stenosis)
- Do not use on patients with solid stool.
- Do not use for patients who have recently (less than 6 weeks ago) had a rectal anastomosis or an anal or sphincter reconstruction
- Do not use for patients with compromised rectal wall integrity (e.g. ischaemic proctitis, mucous membrane laceration)

- Not suitable for paediatric use
- Do not use in the event of pronounced haemorrhoids with corresponding symptoms
- Do not use for patients with a suspected or diagnosed rectal/anal tumour

6. Warnings & precautions

1. Only fill the air quantity specified for use (85 - 90 ml) into the low pressure balloon, do not use an impermissible air quantity (over-filling) or liquids.
2. Any emergent rectal blood must be tested to ensure that there are no signs of decubital gangrene. If there are signs of decubital gangrene, we recommend removing the catheter.
3. Do not use any mechanical pumps to rinse the system.
4. Do not perform irrigation on patients with compromised intestinal wall integrity.
5. Excessive and prolonged pulling on the catheter may cause erosion to the anal canal or anoderm, or can cause the low pressure balloon to drift in the anal canal.
6. Patients with a weak or too strong a sphincter function may eject the catheter when used normally.
7. Take particular care on patients whose rectum has been changed due to a stricture caused by radiotherapy or whose rectum is suffering from radiation proctitis.
8. The patient must not lie on the drain hose or connectors for the auxiliary lumen.
9. Do not sterilise and do not re-use.
10. **Diarflex® Pro** is not suitable for use for more than 29 consecutive days.
11. Patients with severe tenesmus or patients who feel tenesmus or severe pain after inserting the catheter may not tolerate the catheter. The catheter must be removed in this case.

7. Side effects

The use of a rectal stool drainage may cause the following unwanted side effects:

- Perforation
- Decubital gangrene
- Loss of sphincter tone
- Obstruction
- Infection
- Excess stool escape

8. Use

8.1 Preparation

In order to ensure optimised catheter performance, all shaped stool must be removed from the intestine prior to inserting the catheter. Furthermore, all stool must be removed from the rectum before inserting the catheter. The preferred patient position for inserting the catheter is the left lateral knee-chest position. Other positions are possible depending on the patient's condition. The aim of the position is to achieve maximum sphincter relaxation for easier catheter insertion.

8.2 Inserting the catheter

1. **(fig. 1)** - Attach the catheter hose with screw cap to the collection bag.
2. **(fig. 2)** - Using a Luer syringe, remove all the air from the low-pressure balloon via the connection labelled "AIR" (then remove the syringe).
3. **(fig. 3)** - Coat the low pressure balloon and the rectal opening generously with water-soluble lubricant.
4. **(fig. 4)** - Grip support shaft
5. **(fig. 5)** - Press the support shaft together in the centre
6. **(fig. 6)** - Fold the outer sides upwards
7. **(fig. 7)** - Grasp the catheter directly behind the folded structure.
8. **(fig. 8)** - Insert the catheter into the rectum with a twisting motion. Note the end of the support shaft for the maximum insertion depth. (The 2nd balloon remains outside the rectum)
9. **(fig. 9)** - Fill the low-pressure balloon with 85 - 90 ml of air using a Luer syringe via the connection labelled "AIR" (then remove the syringe).

10. **(fig. 10)** - Ensure positioning is correct

11. **(fig. 11)** - A catheter syringe (not included in the scope of delivery) can be used to take a faecal sample via the sampling port.

1. *Moisten the tip of the catheter syringe*
2. *Open sampling port*
3. *Remove a sufficient amount of stool*
4. *Close the sampling port again*

8.3 Caring for and monitoring the system

- At regular intervals, check that the catheter and the collection bag are positioned so that the catheter hose is not twisted, kinked or pressed together externally.
- Check at regular intervals that the catheter is not migrating into the bowel (the end of the support shaft with balloon should always be visible).
- At regular intervals, check that no stool collects in the catheter's drain hose. (To remove this, you can use your hand to perform a guided pressing movement to push the stool through the discharge hose towards the collection bag from the outside)
- Ensure and check at regular intervals to ensure that the patient is not lying on the drain hose or the connectors for the auxiliary lumen.
- Rinse the inner lumen of the catheter hose at least twice per day (more frequently if necessary) in order to remove the collected stool quantity. To do this, inject water via the connector labelled "Irrigation" using a Luer syringe (then remove the syringe). Allow the rinsing water to run into the collection bag.
- When the collection bag is full, replace it with a new bag. For this, the hose can be temporarily sealed with one of the enclosed sealing plugs. Dispose of the used sealing plug afterwards.
- Close the full collection bag with the attached screw cap and dispose of it.
- To dispose of the hose at the end of the treatment, it is closed with the sealing plug attached to the end of the hose and disposed of.

8.4 Administering medication

It is possible to administer medication rectally via the connector labelled "Irrigation".

Attention!

This may only be done under doctor's orders and the medication dose must be checked by the doctor!

- Visually check whether the irrigation line is free, flush via the connector labelled "Irrigation" if necessary
 - Prepare the Luer syringe with the medication per doctor's orders
 - Place the enclosed catheter clamp approx. 10 cm below the outer balloon, but do not close it
1. **(fig. 12)** - Attach the Luer syringe to the connector labelled "Irrigation" and administer the medication. Then rinse the catheter immediately with 50 ml tap water (body temperature) via the same connector so that the medication enters the rectum.
 2. **(fig. 13)** - Close the catheter clamp.
 3. **(fig. 14)** - After a medication exposure time specified by the doctor, re-open the catheter clamp and remove it.

8.5. Removing the catheter

1. **(fig. 15)** - Remove all the air from the low-pressure balloon via the connection labelled "AIR" using a Luer syringe. Then remove the syringe.
2. **(fig. 16)** - Apply water-soluble lubricant to the rectal opening, grasp the catheter just behind the rectal opening and pull the low pressure balloon out carefully and with even pressure.

Attention!

Due to the catheter's elasticity, it may be ejected from the anus very suddenly. Please take care to protect yourself from possible contamination.

9. Disposal

The product may only be disposed of in accordance with the applicable national regulations for waste materials.

10. Reporting

All serious incidents that occur in conjunction with using the product must be reported to the manufacturer and the authority responsible.

11. Warranty

We guarantee that the delivered item is free of defects in accordance with our general terms and conditions (T&C). A copy of our T&C is provided on our business documents and on our website at www.primed-halberstadt.de. A failure to comply with the instructions voids the warranty!

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