

Medical aids for  
colon diagnostics

# PRODUCT CATALOG



**Primed<sup>®</sup>**





**Primed<sup>®</sup>** Medizintechnik GmbH is a highly experienced manufacturer of medical devices.

For many years we have been supplying our customers with high quality products "MADE IN GERMANY".


**Primed<sup>®</sup>** both develops and manufactures all products in-house, in a certified development process according to DIN EN 13485.

Starting from the product idea, through material tests, tool construction and product release, we produce innovative medical devices from one single source. Our production covers the entire manufacturing process (injection moulding/extrusion/assembly), up to sterilization in our own EO sterilization facility.

These conditions enable us to react extremely flexible to market requirements and to ensure the continuous development of our products.

- Rectal catheters.....5
- Rectal catheters without balloon.....6
- Rectal catheters with balloon .....7
- Specialized catheters ..... 10
- Contrast medium bags ..... 12
- Accessories..... 16

Rectal catheters are used in colon diagnostics during radiological examinations. They are applied in mono or double contrast X-ray procedures, enabling a thorough examination of the colon. Intestinal tubes are available according to patient's requirements, with or without balloon and with or without integrated air insufflation tube.



**RECTAL CATHETERS WITHOUT BALLOON**



**Rectal catheter for adults**

- one-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- central orifice and 4 lateral eyes
- with transparent stepped connector
- length = 160 mm

PU/ SU

ART.NO.

10/100

20073



**Rectal catheter for adults,  
with air supply (double contrast)**

- one-piece, with olive-shaped tip (Ø 17 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- central orifice and 2 lateral eyes
- tube connection for air insufflation
- with transparent stepped connector
- length = 180 mm

PU/ SU

ART.NO.

10/100

21548

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

**RECTAL CATHETERS WITH BALLOON**



**Rectal balloon catheter for adults**

- one-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with Luer-Lock connector
- central orifice and 2 lateral eyes
- with transparent stepped connector
- including Luer-Lock syringe
- length = 155 mm

PU/SU

ART.NO.

10/100

21546



**Rectal balloon catheter for adults**

- one-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with Luer-Lock connector
- central orifice and 2 lateral eyes
- with transparent stepped connector
- length = 155 mm

PU/SU

ART.NO.

250/250

20101

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

**RECTAL CATHETERS WITH BALLOON**



**Rectal balloon catheter for adults**

- two-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with Luer-Lock connector
- central orifice and 2 lateral eyes
- with transparent stepped connector
- including Luer-Lock syringe
- total length = 380 mm (including tube extension)

PU/SU

ART.NO.

10/100

21545



**Rectal balloon catheter for adults**

- two-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with Luer-Lock connector
- central orifice and 2 lateral eyes
- with transparent stepped connector
- total length = 380 mm (including tube extension)

PU/SU

ART.NO.

125/125

20089

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT



**RECTAL CATHETERS WITH BALLOON**



**Rectal balloon catheter for adults**

- two-piece, with olive-shaped tip (Ø 16 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with funnel connector
- central orifice and 2 lateral eyes
- with transparent stepped connector
- total length = 380 mm (including tube extension)

PU/SU

ART.NO.

125/125

20133



**Rectal balloon catheter for adults, with air supply (double contrast)**

- one-piece, with olive-shaped tip (Ø 15 mm)
- conical shank made of soft, flexible PVC (DEHP free)
- soft silicone balloon with Luer-Lock connector
- central orifice and 2 lateral eyes
- tube connection for air insufflation
- with transparent stepped connector
- including Luer-Lock syringe
- length = 180 mm

PU/SU

ART.NO.

10/100

21544

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

SPEZIALIZED CATHETERS



**Rectal catheter for children**

- one-piece, without olive-shaped tip
- conical shank made of soft, flexible PVC (DEHP free)
- central orifice and 2 lateral eyes
- atraumatic tip (Ø 8 mm)
- with transparent stepped connector
- length = 170 mm

PU/SU

ART.NO.

10/20

21547



**PRIREKTOSILK, rectal catheter for neonates**

- conically shaped anal olive with handle (length = 3.5 cm)
- supply with Luer-Lock connection (length = 30 cm)
- a specially shaped olive enables secure sealing of the anal orifice
- made of soft, flexible medical-grade silicone

PU/SU

ART.NO.

10/100

20693

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

SPECIALIZED CATHETERS



**Anus Praeter catheter**

- made of soft, flexible PVC (DEHP free)
- for radiological examinations of patients with artificial bowel outlet
- tube catheter with lateral orifice at the front end and stepped connector at the back end
- with conical bulge for secure sealing of the anus and adhesive plaster for fixing on the skin
- with slide clamp

PU/SU


ART.NO.

1/125

21543

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

We provide a comprehensive range of contrast medium bags.



**CONTRAST MEDIUM BAGS**



**Contrast medium bag**

- 2500 ml bag with sealing cap
- firmly attached tube with one-handed clamp (length = 180 cm)

PU/ SU

ART.NO.

80/80

21424



**Contrast medium bag  
with stop ball**

- 2500 ml bag with stop ball and sealing cap
- firmly attached tube with one-handed clamp (length = 180 cm)

PU/ SU

ART.NO.

80/80

21425

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

**CONTRAST MEDIUM BAGS**



**Contrast medium bags  
with air supply**

- 2500 ml bag with Luer-Lock valve and sealing cap
- firmly attached tube with one-handed clamp (length = 180 cm)

PU/SU

ART.NO.

80/80

21421



**Contrast medium bag  
with Y-piece**

- 2500 ml bag with sealing cap
- firmly attached tube with one-handed clamp and Y-piece (length = 200 cm)
- Y-piece enables attachment of additional bag and can be closed individually with one-handed clamp

PU/SU

ART.NO.

80/80

21422

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

CONTRAST MEDIUM BAGS



Contrast medium bag  
with screw cap

- 2000 ml bag with screw cap
- firmly attached tube with slide clamp (length = 150 cm)

PU/SU

ART.NO.

80/80

21541

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT

## PUMP &amp; SEAT PAD

Pump for air insufflation

- for air insufflation into intestine during double contrast procedure
- with black stepped connector for attaching blue tube of double contrast rectal catheters with air supply (item codes 21548, 21544)
- latex-free
- DEHP-free

PU/SU

ART.NO.

10/300

21549

Seat pad for contrast enema

- inner Ø 34.0 cm
- easily inflatable
- latex-free
- DEHP-free

PU/SU

ART.NO.

1/25

21557

PU/ SU = PACKAGING UNIT/ SHIPPING UNIT





MANUFACTURER

---

**Primed**<sup>®</sup> HALBERSTADT  
MEDIZINTECHNIK GMBH  
STRASSE DES 20. JULI 1  
D-38820 HALBERSTADT  
PHONE: 0049-(0)3941-668-6  
FAX: 0049-(0)3941-24565  
PRIMED@PRIMED-HALBERSTADT.DE

DISTRIBUTOR

---

ALL RIGHTS TO TECHNICAL MODIFICATIONS RESERVED!

REVISION STATUS: JULY 2015